

RheoSense microVISC SOP

Equipment

RheoSense μ VISC



Contacts

Graduate Student: Katie Copenhaver (katie.copenhaver@gtri.gatech.edu)

Russo Lab Associate: Alyssa Blake (ablake8@gatech.edu)

Undergraduate Student Contact: Ian Watt (iwatt3@gatech.edu)

Why Should I use this Equipment?

To quickly, easily, and accurately determine viscosities of solutions at or near room temperature using a small sample volume.

Safety tips for user and instrument

- 1) Prepare your sample far enough away from the machine to avoid chemical spills, especially on connectors.
- 2) Always use gloves when handling samples or system components.
- 3) After drawing a sample into the pipette, wipe the outside dry with a lint-free wipe before inserting into the microVISC.
- 4) Never manually depress a pipette plunger after you mount it into the unit (could result in permanent damage to the sensor).
- 5) Do not allow particles larger than 10% of the sensor channel depth in the sample. A table of chip types and operating ranges is available in the appendix. Check the chip type before taking measurements to ensure that your sample can be safely and accurately measured.
- 6) Always perform a cleaning cycle with appropriate solution before and after measurements.
- 7) Samples above 100cP will take more than one cleaning cycle.

- 8) When using acidic or basic solvents, try to take measurements without breaks between runs and run a cleaning cycle immediately after measurements are finished. Prolonged exposure to acids or bases can damage the sensor.

Quick start user guide

- 1) Enable temperature control using on/off button on the left side of the panel.
 - a) Change the set point using the arrows (tap for 0.1° increments, hold for 1° increments) or by using the numeric keys when the set point is flashing.
 - b) Advanced options (temperature margin, soak time, etc) can be changed by pressing “Menu” and “2” for temperature control.
- 2) Load a pipette with sample by grasping the flange and slowly drawing back the plunger while tip of pipette is in the sample
- 3) Mount the pipette in the unit by guiding the tip into the sample inlet port and firmly pushing down on the flange until it clicks into place
- 4) Press the RUN/STOP button on the keypad to run in Auto, Advanced, or Cleaning mode. The mode can be changed by pressing ENTER on the main panel and using the arrow keys.
- 5) Remember these three rules
 - a) MNF (measure nothing first): in this case, “nothing” is a 1% Aqueous solution in the cleaning cycle. A cleaning cycle **must** always be run before **and** after sets of measurements.
 - b) MSS (measure something simple): in the case of polymer solutions or suspensions, measure the pure solvent.
 - c) MYS (measure your stuff): if measuring samples at a range of concentrations/viscosities, start with the lowest and work your way up.

Helpful hints

- 1) Check the sensor type before beginning a measurement to ensure that the viscosity of your sample will be in range of the sensor’s measurement capacity.
- 2) Always hold pipettes by the rectangular flange to minimize heat transfer to the sample.
- 3) Always try to load the pipette so there are no bubbles in it.
 - a) After loading the pipette with sample, hold the pipette vertically and push the plunger up slightly until the bubble is expelled from the tip.
- 4) If possible, load your sample before enabling the temperature control. After a new sample is loaded, you must wait for the temperature to re-stabilize, which can take 3-5 minutes.
- 5) If repeated measurements of one sample show a gradual increase or decrease in viscosity, the priming volume may have been insufficient. Repeat the measurement until the reading is stable.
- 6) To conserve sample volume when measuring consecutive samples with large differences in viscosity, run a cleaning cycle between measurements.

Operating Modes

- 1) Automatic
 - a) After loading the sample and pressing RUN/STOP, the system will automatically move the pusher to the plunger position, prime the sensor with test fluid (10-200µL), pause for

sample relaxation (10-20s), take a baseline measurement, and measure the viscosity of your sample.


- b) During the measurement, the following values are displayed:
 - i) Viscosity (mPa-s)
 - ii) Temperature (°C)
 - iii) Pumping rate (s⁻¹ or μL/min)
 - iv) Dispensed volume (μL)
 - v) Pressure (%age of full scale)
 - vi) Measurement progress bar
 - c) Other measurement parameters can be viewed using the up/down arrow keys.
 - d) After a measurement is complete there are 3 options:
 - i) Press RUN/STOP to take another measurement of the same sample.
 - ii) Press MENU to access advanced parameters and functions.
 - iii) Press .HOME to return pusher to home position and enable removal of the sample.
- 2) Advanced
- a) In advanced mode, you can specify all of the following parameters (or leave some on auto for the measurement system to decide):
 - i) Rate: units are either in s⁻¹ for shear rate or μL/min for flow rate. Flow rate range is 0.5-450 μL/min.
 - ii) Measurement volume: available measurement volume is 5-400 μL.
 - iii) Priming volume: available priming volume is 1-200 μL.
 - iv) Pause time (for stress dissipation after priming): available pause time is 0-36000 seconds.
- 3) Cleaning
- a) To clean the sensor cartridge, fill a new pipette with solvent that is miscible with the samples used. The system will dispense the entire contents of the pipette at 100 μL/min. After running a cleaning cycle of solvent, repeat with 1% Aqueous solution.

Appendix


1) Chip types and operating ranges

Chip type	Viscosity range (cP)	Flow channel depth (μm)
A05	0-100	50
A10	4-600	100
A20	10-2,000	200
A30	30-8,000	300
B10	60-5,000	100
B20	70-10,000	200
B30	100-20,000	300
C10	400-22,000	100
C20	500-40,000	200
C30	2,000-80,000	300


2) Hardware Overview



3) Main Panel Overview



4) Menu Structure



5) LCD Display Overview

